

EA-100

DEMONSTRATION

“QUICK START”

CASIO DATA COLLECTOR-“Manual Set-up”

NO PROGRAMMING NECESSARY! NO CALCULATOR NECESSARY

**PURPOSE: TO MEASURE THE CHANGE IN TEMPERATURE OVER A SPECIFIED TIME
(Using your hand)**

STEPS FOR “MANUAL-MODE” (No program)

Make sure temperature probe is connected to EA-100. Do not connect to calculator.

You will program it by choosing 3 built in options:

	<u>OPTION RANGE/CHOICES</u>
1) Sample Time	10 milliseconds to 60 seconds
2) Number of samples to be collected	10 to 200
3) Time consideration...	(0) none, (1) absolute, or (2) relative

<u>KEYSTROKE</u>	<u>EXPLANATION</u>
1) Press [SHIFT] [MODE]	“Set-up”
2) Press [DataLOG] until “1.00sec” appears see the	“Next” Press [DataLOG] several times until you sample time you want to use on your screen. In this case 1.00sec.
3) Press [TRIGGER]	Chooses 1.00sec as sample time.
4) Press [DataLOG] until “10” appears	Press [DataLOG] several times until you see the number of samples you want to use on your screen. In this case 10.
5) Press [TRIGGER]	Chooses 10 as number of samples to collect.
6) Press [DataLOG] until “1” appears	“Next” Press [DataLOG] until you see the time consideration you want to use on your screen. In this case 1. 1 means time is collected as data and in this case time is stored as,: 1, 2, 3, 4, 5, 6, 7, 8, 9, 10.
7) Press [TRIGGER]	Chooses 1 as time consideration (absolute time).

ARE YOU READY???????

Put the temperature probe in your hand.

8) Press [TRIGGER] Starts data collection.

WAIT UNTIL SCREEN SAYS “DONE”!

*Now view your data by pressing [DataLOG] key over and over again.

*To redo or repeat the last experiment....press [SHIFT] [RESTART].

*To transfer the data from the EA-100 to calculator...Use this program (or type this in the RUN icon)...

Receive(List1) [EXE] Receive (List 2) [EXE]

Data will show up in list 1 and list 2 in statistics icon.

*You may use any other probe with this manual set-up, not just the light probe.

PATH OF A BALL

MATERIALS:

- 1) EA-100
- 2) MOTION DETECTOR
- 3) GRAPHING CALCULATOR WITH SONIC PROGRAM
- 4) 1 LARGE "BOUNCY" BALL

PROCEDURE:

Measure the track of a ball dropped, based on it's distance from the floor

- 1) Connect the motion detector to SONIC CH.
- 2) Connect the data collector to the CFX calculator. Make sure it's plugged all the way in.
- 3) Turn the data collector on and execute the program "Sonic" on the CFX.
- 4) Secure the motion detector above where the ball will be dropped and get the ball ready to drop under the detector.
- 5) Push "trigger" on the data collector at the same time as dropping the ball to start the experiment.
- 6) Measurement will be finished when "done" is displayed on the data collector.
- 7) Push "EXE" on the CFX to graph the results.

NOTES:

-Alter the sample interval if you like...

Enter the program icon on the CFX-9850 and highlight the "sonic" program.

[F2] to edit program.

GO to 6th line and alter it as you wish following the example below.

Example {3, .02, 150, 1}

{3 means motion detector, .02 is time interval, 150 is number of data points to collect, 1 is absolute time}

QUESTIONS:

- 1) Predict the graph
- 2) Copy the graph to a piece of paper. Are there any points that don't seem to belong? Circle them.
- 3) What are some reasons for these points?
- 4) How do these points effect the data?
- 5) Describe the path of the ball.
- 6) Is it the same at the beginning and at the end of the data collected?
- 7) What would the graph have looked like if the ball you dropped was bigger?
- 8) What would the graph have looked like if the ball you dropped bounced better?
- 9) How was the actual graph different from the one you predicted? Why was it different?

WHICH LIGHT SHOULD YOU BUY?

MATERIALS:

- 1) EA-100
- 2) LIGHT PROBE
- 3) 2 different brand flashlights or light bulbs
- 4) measurement stick
- 5) masking tape

PROCEDURE:

Measure the light given off by each flash light at distances of every 6 inches.

- 1) Connect the light probe to CH 1.
- 2) Turn the data collector on.
- 3) Use [shift] [setup] to choose 5 second intervals, and 10 data points, and 1 for time.
- 4) Place 1 flash light on a table and use the masking tape to make a mark every 3 inches for 2 1/4 feet.
- 5) Turn the flash light on and put the light detector at the 0" mark. (Is it in the center of the light beam?)
- 6) Push "trigger" on the data to start the experiment.

MOVE the light probe to the next mark after each data point is collected. You have a 5 second interval to move to the next mark before the next point will be collected.

You should finish on the 2' 3" mark when EA-100 reads "done".

NOTES:

HI - means reading too high for probe to measure...out of range.

LO - means reading too low for probe to measure...out of range.

Secure light probe in a place where the "center" of the light beam is directly across from probe.

QUESTIONS:

- 1) Predict which light is better? Why?
- 2) Describe different ways that would make one light better than another.
- 3) Copy the data points to the table below.

REPEAT THE EXPERIMENT FOR FLASHLIGHT #2. Copy the data points to the table.

- 1) Both flash lights were sold for \$5 at the local store. Which one is a better buy (gives off the most light)?
- 2) Describe why you made your decision.

	LIGHT # 1	LIGHT # 2
0"		
3"		
6"		
9"		
12"		
1'		
1' 3"		
1' 6"		
1' 9"		
2'		

HEATING AND LINE OF BEST FIT

MATERIALS:

- 1) EA-100
- 2) TEMP PROBE
- 3) YOUR HAND

PROCEDURE:

Measure the heating of a temperature probe in your hand, transfer data to the calculator, and find the line of best fit.

- 1) Connect the temperature probe to CH 1.
- 2) Turn the data collector on.
- 3) Use [shift] [setup] to choose 1 second intervals, and 10 data points, and 1 for time.
- 4) Push "trigger" on the data to start the experiment when you put the probe into your hand.
- 5) EA-100 will say "done" when complete.
- 6) Connect EA-100 to calculator and execute "RECEIVE1" program on the CFX-9850.
- 7) Make sure the cable is plugged all the way into both units. Press until you hear a click.
- 8) Go to the Stat Icon on the calculator and graph your results.
- 9) While viewing your graph choose a line of best fit.

QUESTIONS:

- 1) Predict the graph you will see.
- 2) Describe your graph.
- 3) Which of the lines of best fit (regressions) is best? (Try each of the F1-F5) that will appear on the bottom of your graph after you draw it. Hint: which one has the highest "r" value?
- 4) How was the actual graph different from the one you predicted? Why was it different?

REAL TIME...CAN YOU MATCH THESE CURVES?

MATERIALS:

- 1) EA-100
- 2) MOTION DETECTOR
- 3) GRAPHING CALCULATOR WITH "REALTIME" PROGRAM

PROCEDURE:

Measure the movement in front of the motion detector in real time. The data collector (EA-100) can measure the variation and graph the results as they are received on the CFX-9850.

- 1) Connect the motion detector to SONIC CH.
- 2) Connect the data collector to the CFX calculator. Make sure it's plugged all the way in. Push hard into both units.
- 3) Put the motion detector on a table in front of you and get ready to move in front of the motion detector.
- 4) Turn the data collector on and execute the program "REALTIME" on the CFX.
- 5) Measurement will be finished when "done" is displayed on the data collector.

NOTES:

Use your body to make curves.

You can press [EXE] to restart program, over and over again.

QUESTIONS:

- 1) Predict how you will draw each graph using your body and write it in the space below.
- 2) Draw the following graphs with your movement in front of the motion detector.
- 3) Describe the movement needed to draw each graph.
- 4) How did you make the actual graph? Was your actions different from the one you predicted? Why was it different?
- 5) How would you draw a vertical line.

Horizontal line

downward sloping line

upward sloping line

parabola opening up

parabola opening down

sin wave curve

COOLING OF A LIQUID

MATERIALS:

- 1) EA-100
- 2) TEMPERATURE PROBE
- 3) GRAPHING CALCULATOR WITH "TEMP" PROGRAM
- 4) Hot cup of coffee, water, or other liquid
- 5) ice cubes

PROCEDURE:

Measure the cooling of a cup of coffee when an ice cube is added.

- 1) Connect the temperature probe to CH 1.
- 2) Connect the data collector to the CFX calculator. Make sure it's plugged all the way in. Push hard into both units.
- 3) Turn the data collector on and execute the program "temp" on the CFX.
- 4) Put the temperature probe into the coffee and wait for the temperature to peak or level out.
- 5) Push "trigger" on the data to start the experiment and then, immediately add 1 ice cube.
- 6) Measurement will be finished when "done" is displayed on the data collector.
- 7) Push "EXE" on the CFX to graph the results.

NOTES:

If you don't see your data on your graph, make sure your "setup" has "AUTO" selected for "Stat Graph".

Do this by go into the STAT icon and press... [SHIFT], [MENU], [F1], [EXIT]

Now re-draw your stat graph 1 using [F1: GRPH 1].

QUESTIONS:

Copy the graph to a piece of paper.

- 1) Predict the graph you will see.
- 2) Describe your results.
- 3) Does the cooling occur at the same rate through out your experiment?
- 4) How was the actual graph different from the one you predicted? Why was it different?

REPEAT THE EXPERIMENT WITH A NEW CUP OF COFFEE

THIS TIME, IN STEP #5, ADD 4 ICE CUBES

- 1) How did your graph change?
- 2) Describe why the graph changed.
- 3) How would the graph change if the sample interval was changed to every 5 seconds, over 10 minutes?
- 4) How would the graph change if the sample interval was changed to every 5 milliseconds, over 2 seconds?